

UPLSAPO Series

Thanks to the application of Olympus' original UW multicoating, these Super Apochromat objectives fully compensate for both spherical and chromatic aberrations from the UV to the near infrared region. Their sensitivity to fluorescence emissions ensures the acquisition of sharp, clear images, without color shift, even in bright field and Nomarski DIC observations. For quality and performance, they offer an unbeatable solution to every kind of digital imaging need.



[UPLSAPO 4X](#)



[UPLSAPO 10X2](#)



[UPLSAPO 20X](#)



[UPLSAPO 20XO](#)



[UPLSAPO 30XS](#)





UPLSAPO 40X2



UPLSAPO 60XO



UPLSAPO 60XW



UPLSAPO 60XS



UPLSAPO 100XO

